SODIO

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Enzymatic Method

Liquid Reagents ready to use

REF. 0237 2x60 ml + 2 ml Calibrators
REF. 0237/2 1x60 ml + 1 ml Calibrators

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE


Quantitative determination of sodium in serum.

CLINICAL SIGNIFICANCE

Measurements obtained by this device are used in the diagnosis and treatment of aldosteronism (excessive secretion of the hormone aldosterone), diabetes insipidus (chronic excretion of large amounts of dilute urine, accompanied by extreme thirst), adrenal hypertension, Addison's disease (caused by destruction of the adrenal glands), dehydration, inappropriate antidiuretic hormone secretion, or other diseases involving electrolyte imbalance.

PRINCIPLE

Sodium is determined enzymatically via sodium-dependent β -galactosidase activity with ONPG as the substrate. The absorbance at 405 nm of the product O-nitrophenyl is to the sodium concentration.

ONPG = o-nitrophenyl- β -D-galactopyranoside

SAMPLE

Serum.

KIT COMPONENTS

Reagent (A) Sodium Volume = 40 ml	Good's Buffer - pH 8.5 Encrypting (> 0.4 mM), β -D-galactosidase (< 8 U/mL), Proclin 300 (0.02%)
Reagent (B) Sodium Volume = 20 ml	Good's Buffer - pH 6.5 O-nitrophenyl β -D-glycoside (> 0.5 mM), Proclin 300 (0.02%)
Low Calibrator Volume = 1 or 2 ml	Buffered Sodium (batch-specific value indicated on the bottle)
High Calibrator Volume = 1 or 2 ml	Buffered Sodium (batch-specific value indicated on the bottle)

Reagents and calibrators are stable until the expiration date indicated on the label if stored at 2-8°C. Do not use over the expiry date.

Once opened reagents are stable for 2 months if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagents, ready to use.

PRECAUTIONS AND WARNINGS


Use the normal precautions required in the laboratory.

The Reagents R(A) and R(B) contain Proclin 300. Avoid ingestion or contact with the skin or mucous membranes. After contact with skin, wash immediately with plenty of water. In case of eye contact or ingestion, seek medical advice immediately. The material safety data sheet are available on request.

Dispose of waste according to local laws.

PROCEDURE

Diagram of the test to automated chemistry analyzers:

CALIBRATION

For the calibration is recommended to use the calibrators High and Low. It is advisable to perform a 2-point calibration every week.

This test uses a linear calculation and a blank reagent.

REFERENCE RANGE

136 – 146 mmol/L (313 - 336 mg/dl)

Conversion factor: mmol/l x 2.3 = mg/dl

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Linearity: the method is linear for sodium concentrations between 80 and 180 mmol/L (184 and 414 mg/dl).

Precision intra-assay:

	Level 1	Level 2
Mean (mM)	128.94	155.84
SD	1.57	1.72
CV %	1.2	1.1

Precision inter-assay:

	Level 1	Level 2
Mean (mM)	128.94	155.84
SD	2.01	2.56
CV %	1.56	1.65

Interferences: The following substances, which are normally present in the serum, have given rise to a deviation less than 10% at the indicated concentrations: NH_4Cl 1.5 mM, KPi 2.0 mM, CaCl_2 7.5 mM, KCl 10 mM, CuCl_2 0.5 mM, ZnCl_2 0.5 mM, FeCl_3 0.5 mM, glucose to 5 mM, ascorbate 10 mM, bilirubin to 40 mg/dl, conjugated bilirubin to 40 mg/dl, hemoglobin to 500 mg/dl and triglycerides to 1000 mg/dl.

Correlation against a reference method: $Y = 1.05x - 2.23$ $r = 0.98$

REFERENCES

- Berry, M.N. et al., (1988) Clin. Chem. 34, 2295
- Tietz, N. W. (1983) Clinical guide to Laboratory Tests, p. 384 W.B. Saunders Co., Philadelphia.