TRANSFERRINA (TIBC-LIBC)

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$Y = 0.962x + 10.28 \quad r = 0.976$

BIBLIOGRAFIA

1. Ramsey W.N.M., The determination of the TIBC of serum, Clin. Chem. Acta 2.221 (1957).
2. Jung D.H., Parek A.C., Am.J. Clin. Path. 54, 813 (1970).
3. Duffy J.R., Gaudin J., Clin. Biochem., 10, 122 (1977).
4. Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods, Philadelphia, W.B. Saunders, P. 1434 (1984).
5. Levy A., Vitacce P., Direct determination and binding capacity serum iron, Clin. Chem. 7, 241–248 (1961).
6. Schade et al., Bound iron and unsaturated iron binding capacity of serum, rapid and reliable quantitative determination, Proc. Soc. Exp. Med. 67, 442 (1954).
7. Young D.S., et al., Clin. Chem. 21:1D (1975).

TRANSFERRIN (TIBC-LIBC)

Colorimetric Method with precipitation

REF. 0022 1x100 ml

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE

Quantitative determination of TIBC e LIBC in serum.

PRINCIPLE

Transferrin is a protein able to bind the iron present in serum. To define the quantity of this protein, the serum is saturated with an excess of trivalent iron ions and then this excess is eliminated by precipitation with basic carbonate magnesium. The transferrin-bound iron is dosed on the supernatant. The transferrin-bound iron is dosed on the supernatant. The transferrin-bound iron represents the total capacity to bind iron TIBC (Total Iron Binding Capacity). The LIBC (Latent Iron Binding Capacity) or free transferrin is obtained by deducting the serum iron value from the TIBC.

SAMPLE

Unhemolyzed fresh serum. Remove serum from clot as soon as possible. Iron in serum is stable 4 days at 15-25 °C and one week at 2-8 °C.

KIT COMPONENTS

Reagent (A) TIBC Saturating solution Volume = 100 ml	Iron chloride	50 mmol/l
Reagent (B) TIBC Precipitant Reagent powder	Magnesium carbonate	21 g
Measuring spoon	1	

Optional: Measuring cup (Ref. 0022/1)

The reagents are stable until the expiration date indicated on the label if stored at room temperature (15-25 °C) and protected from light.

Do not freeze. Once opened reagents are stable for 2 months at 15-25°C if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

The reagents are ready to use.

PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Add 1 ml of Reagent (A) to 0.5 ml of sample in centrifuge cuvette.

Mix on vortex. Leave the solution stand for 5 minutes and then add 2 measuring spoons of Reagent (B), leave at room temperature for 10-20 minutes, mixing 3 or 4 times on vortex.

Then centrifuge for 5 minutes at about 3000 rpm/min.

Determine iron concentration in the limpid supernatant.

For this purpose we recommend the use of the kit IRON F (Ref. 0086 – 0089).

RESULTS CALCULATION

$TIBC (\mu\text{g/dl}) = \frac{A_x - A_{bx}}{A_s} \times \text{Standard Value} \times 3$ (serum dilution)

$TIBC (\mu\text{g/dl}) = \text{Iron concentration in the supernatant} \times 3$ (serum dilution)

$LIBC = TIBC - \text{Serum Iron}$

Conversion factor: $\mu\text{g/dl} \times 0.1791 = \mu\text{mol/l}$

EXPECTED VALUES

TIBC: 250 – 420 $\mu\text{g/dl}$ (45 – 75 $\mu\text{mol/l}$)

LIBC: 150 – 340 $\mu\text{g/dl}$ (27 – 61 $\mu\text{mol/l}$)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions.

PERFORMANCE

Sensitivity: the sensitivity of the method is: 5 $\mu\text{g/dl}$.

Linearity: the method is linear up to 1000 $\mu\text{g/dl}$. For higher values, dilute the sample 1:2 and multiply the result by 2.

Precision intra-assay:

	Level 1	Level 2
Mean ($\mu\text{g/dl}$)	279	392
SD	3.06	3.14
CV %	1.3	1.5

Precision inter-assay:

	Level 1	Level 2
Mean ($\mu\text{g/dl}$)	287	390
SD	3.36	5.73
CV %	5.7	4.2

Interferences: bilirubin does not interfere up to 55 mg/dl. Hemoglobin does not interfere up to 200 mg/dl. Copper does not interfere up to 500 $\mu\text{g/dl}$.

Correlation against a reference method:

$Y = 0.962x + 10.28$ $r = 0.976$

REFERENCES

1. Ramsey W.N.M., The determination of the TIBC of serum, Clin. Chem. Acta 2.221 (1957).
2. Jung D.H., Parek A.C., Am.J. Clin. Path. 54, 813 (1970).
3. Duffy J.R., Gaudin J., Clin. Biochem., 10, 122 (1977).
4. Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods, Philadelphia, W.B. Saunders, P. 1434 (1984).
5. Levy A., Vitacce P., Direct determination and binding capacity serum iron, Clin. Chem. 7, 241–248 (1961).
6. Schade et al., Bound iron and unsaturated iron binding capacity of serum, rapid and reliable quantitative determination, Proc. Soc. Exp. Med. 67, 442 (1954).
7. Young D.S., et al., Clin. Chem. 21:1D (1975).