FERRO F

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IRON F

Colorimetric Method with liquid Ferene

REF. 0089 4x50 ml
REF. 0086 8x50 ml

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE

Quantitative determination of iron in the serum.

CLINICAL SIGNIFICANCE

In biological fluids, iron is bound to proteins such as hemoglobin, myoglobin, transferrin, and ferritin. Elevated iron levels may indicate conditions such as hemochromatosis and liver damage, while low levels may be a sign of iron-deficiency anemia due to malabsorption related to gastrointestinal disorders or heavy menstrual bleeding. For a more detailed overview of iron metabolism, it is also recommended to determine transferrin and ferritin levels.

PRINCIPLE

Serum iron bound to transferrin is released in acid environment. Fe (III) ions are then reduced to Fe (II), which binds to the Ferene to give a stable colored complex, whose color intensity is proportional to the amount of iron in the sample.

SAMPLE

Fresh serum. Do not use hemolyzed samples.

Separate serum from clot as soon as possible.

The iron in the sample is stable 4 days at room temperature and at least one week at 2-8°C.

KIT COMPONENTS

Reagent (A) Fe F Volume = 45 ml	Buffer Complexing Guanidine hydrochloride Tensioactive agents	400 mmol/l > 60 mmol/l 1 M	
Reagent (B) Fe F Volume = 20/40 ml	Ferene Reducing Thiourea	31 mmol/l 17 mmol/l > 200 mmol/l > 117 mmol/l	
Standard Fe Volume = 10 ml	Iron	Value on label	

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C. Do not use over expiry date.

Once opened reagents are stable for 2 months at 2-8°C if contamination is avoided.

REAGENT PREPARATION

Bring all reagents at room temperature before use.

Reagents (A), (B) and the Standard are liquid, ready to use.

PRECAUTIONS AND WARNINGS

Reagent (A) – Contains Acetic Acid (CAS 64-19-7)

Signalwords: **Warning**

H319 Causes serious eye irritation.

H315 Causes skin irritation.

P280 Wear protective gloves / eye protection / face protection.

P337+P313 If eye irritation persists: Get medical advice / attention.

P264 Wash with water thoroughly after handling.

Reagent (B) – Contains Thiourea (CAS 62-56-6)

EUH031 Contact with acids liberates toxic gas.

EUH210 Safety data sheet available on request.

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow. Use the normal precautions required in the laboratory. Dispose of waste according to local laws.

PROCEDURE

Wavelength:	593 nm (578 – 700)
Lightpath :	1 cm
Temperature:	37°C
Reading:	against blank
Method:	Blank sample

Reagent use in manual procedure:

pipette:	Blank Reagent	sample	standard
Reagent (A)	900 µl	900 µl	900 µl
water	200 µl		
sample		200 µl	
standard			200 µl

Mix, read the absorbance of blank sample (Abx) against blank reagent. Add in the same tubes:

Reagent (B)	100 µl	100 µl	100 µl
--------------------	--------	--------	--------

Mix, incubate at 37°C for 5 minutes and read the absorbance of the sample (Ax) and the standard (As) against blank reagent.

Reaction volumes can be proportionally varied.

In the manual procedure prepare a series of test-tubes for blanks reading.

For automated procedure, ask for specific applications.

RESULTS CALCULATION

Iron µg/dl = (Ax – Abx)/(As) x (Standard value)

Conversion factor: µg/dl x 0.1791 = µmol/l

EXPECTED VALUES

Men: 59 – 158 µg/dl (10.6 – 28.3 µmol/l)

Women: 37 – 145 µg/dl (6.60 – 26.0 µmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 5 µg/dl.

Linearity: the method is linear up to 800 µg/dl. For higher values, dilute the sample 1:2 and multiply the result by 2.

Precision intra-assay:

	Level 1	Level 2	Level 3
Mean (µg/dl)	70.45	115.7	326.5
DS	1.29	1.16	1.84
CV %	1.83	0.75	0.56

Precision inter-assay:

	Level 1	Level 2	Level 3
Mean (µg/dl)	76.07	167.5	336.0
DS	1.38	1.72	2.45
CV %	1.81	1.02	0.73

Interferences: the copper does not interfere up to 400 µg/dl.

Bilirubin up to 10 mg/dl does not interfere.

Hemoglobin interferes with the test. Do not use hemolyzed samples.

Correlation against a reference method: Y = 1.0099x + 0.982

r = 0.9872

REFERENCES

- Higgins T., Clin. Chem. 27, 1619 (1981).
- Vassault, A. et al. Ann. Biol. Clin., 44,686 (1986).
- Young D. S., et al, Clin. Chem. 21:1D (1975).