CLORURI	
Metodo colorimetrico Mercurio Tiocianato	
Reagente liquido pronto all'uso	
REF. 0020	2x100 ml
REF. 0049	4x100 ml
CE	
IVD	

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CHLORIDE

Mercury Thiocyanate Colorimetric Method
Liquid Reagent ready to use

REF. 0020 2x100 ml
REF. 0049 4x100 ml

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE

Quantitative determination of Chlorides in serum, plasma and urine according to mercury thiocyanate reaction.

PRINCIPLE

The chloride ions present in the sample react with mercuric ions, releasing the same quantity of thiocyanate ions which, in the presence of ferric ions, forms a red colored complex. The color intensity is proportional to the chloride concentration present in the sample.

SAMPLE

Serum, heparinized plasma, urine 24/h diluted 1:2 with distilled water and acidified with 2-3 drops of HCl 23 %.

Chlorides in the sample are stable up to 7 days at 2-8°C and one month at -20°C.

KIT COMPONENTS

Reagent (A) Cl Volume = 100 ml	Mercury thiocyanate Mercury nitrate Iron nitrate Chloridric acid	25 mmol/l 0.2 mmol/l 25 mmol/l 60 mmol/l
Standard Volume = 10 ml	Sodium Chloride Sodium azide	Value on label

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C and protected from light. Do not use over expiry date.

Once opened reagents are stable for 2 months if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagent ready to use.

PRECAUTIONS AND WARNINGS

EC Regulation 1272/2008 (CLP):

Reagent (A):

EUH071 Corrosive to the respiratory tract.

Reagents may contain some non-reactive and preservative components.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength: 480 nm (450 – 500)
Lightpath: 1 cm
Temperature: 25, 30, 37°C
Reading: against blank reagent
Method: Increasing End Point
Sample/Reagent: 1/100

pipette:	blank	sample	standard
Reagent (A)	3000 µl	3000 µl	3000 µl
water	20 µl		
sample		20 µl	
standard			20 µl

Mix, incubate at 25, 30, 37°C for 5 minutes. Read the absorbance of the sample (Ax) and the standard (As) against the blank reagent.

The reaction volumes can be proportionally modified.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific application.

RESULTS CALCULATION

Serum, plasma:

Chloride mEq/l = Ax/As x (Standard Value)

Urine 24/h:

Chloride mEq/l = Ax/As x (Standard Value) x 2 (dilution factor) x Urine Vol 24/h (L)

Conversion factors: mEq/L = mmol/l

mmol/l x 3.545 = mg/dl

EXPECTED VALUES

Serum/plasma:

95 – 110 mEq/L (337 - 390 mg/dl)

Urine:

160 – 250 mEq/24h (5.67 - 8.86 mg/24h)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 1 mEq/L.

Linearity: the method is linear up to 140 mEq/L. For higher values, dilute the sample 1:2 with saline and multiply the result by 2.

Precision intra-assay:

	Level 1	Level 2	Level 3
Mean (mEq/l)	46.89	88.80	117.5
SD	0.25	0.29	0.53
CV %	0.53	0.33	0.45

Precision inter-assay:

	Level 1	Level 2	Level 3
Mean (mEq/l)	46.15	90.54	121.1
SD	0.22	0.55	0.74
CV %	0.48	0.61	0.61

Interferences: bilirubin does not interfere up to 20 mg/dl. Triglycerides do not interfere up to 300 mg/dl. For particularly icteric, hemolyzed and lipemic sample, prepare a blank sample by adding 2000 µl of saline solution to 20 µl of sample. Read at 480 nm the absorbance of the blank sample against a saline solution and then deduct it from the sample absorbance (Ax).

Correlation against a reference method:

$Y = 0.9898x + 1.0063$ $r = 0.9871$

REFERENCES

1. Kaplan LA, Pesce AJ, Clinical Chemistry, Mosby (1989)
2. Schales O. and Shales S., J. Biol. Chem. 140, 879 (1941)
3. Young D.S. et al., Clin. Chem. 21:1D (1975)