ZINCO SL

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ZINC SL

Direct Colorimetric Method (5-Br-PAPS)
Liquid Reagents ready to use

REF. 0233 2x50 ml
REF. 0234 4x50 ml

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE

Quantitative determination of zinc in serum, plasma, urine and semen.

PRINCIPLE

In the presence of 5-Br-PAPS, zinc forms a colored compound whose color intensity is proportional to the zinc concentration present in the sample.

Any interference due to copper, iron, cobalt and other trace elements present in the sample is eliminated using specific masking agents.

SAMPLE

Serum or unhemolyzed plasma. Use only heparin as anticoagulants.

24 hours Urine. Dilute urine samples 1:10 with saline.

Semen. Centrifuge semen for 10 -15 minutes at 3000 rpm. Dilute the supernatant 1:100 with saline.

Always consider the dilution factor in calculations.

KIT COMPONENTS

Reagent (A) Zn Volume = 40 ml	Buffer Surfactant complexing	200 mmol/l
Reagent (B) Zn Volume = 10 ml	Buffer 5-Br-PAPS	200 mmol/l 0.2 mmol/l
Standard Zn Volume = 5 ml	Zinc derivative	Value on label

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C. Once opened reagents are stable for 6 months at 2-8°C if contamination is avoided. Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagents, bring to room temperature before use.

PRECAUTIONS AND WARNINGS

Reagents may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength:	546-565 nm
Lightpath :	1 cm
Temperature:	25/30/37°C
Reading:	against blank reagent
Method:	End Point
Sample/Reagent	1/20

Bireagent Procedure:

pipette:	blank	sample	standard
Reagent (A)	960 µl	960 µl	960 µl
water	40 µl		
sample		40 µl	
standard			40 µl

Mix, incubate for 1 minute at room temperature (20-25 °C), read the absorbance of the sample (Ax1) and the Standard (As1) against blank reagent. Then add:

Reagent (B)	240 µl	240 µl	240 µl
-------------	--------	--------	--------

Mix, incubate for 3 minutes at room temperature (20-25 °C), read the absorbance of the sample (Ax2) and the Standard (As2) against blank reagent.

It is possible the use as monoreagent with limpid sera, maintaining a ratio 4+1, ie 4 parts of Reagent (A) + 1 part of Reagent (B).

If necessary, it is possible to proportionally variate the indicated reaction volumes.

Readings are not possible at 578 nm and at wavelengths other than those recommended.

Reaction volumes can be proportionally varied.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

RESULTS CALCULATION

Serum, plasma:

Zinc (µg/dl) = $Ax2 \cdot Ax1 / As2 \cdot As1 \times$ (standard Value)

Urine:

Zinc (µg/dl) = $Ax2 \cdot Ax1 / As2 \cdot As1 \times$ (standard Value) $\times 10$ (dilution) $\times 24$ h Urine Vol (dl)

Semen:

Zinc (µg/dl) = $Ax2 \cdot Ax1 / As2 \cdot As1 \times 200 \times 100$ (dilution)

Conversion factor: µg/dl $\times 0.153 =$ µmol/l

EXPECTED VALUES

Serum: 70 – 150 µg/dl (10.7 – 22.9 µmol/l)

Urine: 150 – 1200 µg/24h (2.3 – 18.4 µmol/24h)

Semen: 2 – 10 mg/dl (0.3 – 1.5 mmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 5 µg/dl.

Linearity: the method is linear up to 1000 µg/dl. For higher values, dilute with saline 1:2 the sample and multiply the result by 2.

Precision intra-assay:

	Level 1	Level 2
Mean (µg/dl)	95.2	135.7
SD	1.03	3.47
CV %	1.10	2.60

Precision inter-assay:

	Level 1	Level 2
Mean (µg/dl)	94.2	133.4
SD	3.49	3.45
CV %	3.70	2.60

Interferences: the copper does not interfere up to 500 µg/dl. Triglycerides do not interfere up to 2000 mg/dl.

Hemoglobin interfere with the test. do not use hemolyzed sample.

Correlation against a reference method: $Y = 0.902x + 8.81$ $r = 0.966$

REFERENCES

1. Akita Abe, Sumiko Yamashita, Clin. Chem. 35/4, 552 – 554 (1989)
2. Tietz Textbook of Clinical Chemistry 2th Ed, Burtis-Ashwood (1994).
3. R. Homsher and B. Zak, Clin. Chem., 31/8, 1310-1313 (1985).
4. T. Makino, M. Saito, Clin. Chim. Acta, 120, 127-135 (1982).