Catena Leggera Lambda

Turbidimetria

REF. 6714 - 25 +5 ml

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Turbidimetry

REF. 6714 - 25+5 ml

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE

Quantitative determination of Lambda Light Chain in human serum by turbidimetric assay.

PRINCIPLE

The addition, in a precise ratio, of the serum containing the Antigen (Lambda Light Chain) to a solution with the specific antiserum (antibody) form an antigen-antibody reaction which determines a turbidity measurable at 340 nm.

SAMPLE

Fresh serum.

The serum may be stored at 2-8°C for 48 hours or, for a longer period, at -20°C.

KIT COMPONENTS

Reagent (A) Buffer Volume = 25 ml	Goods buffer, PEG. Sodium azide 0.95 g/l
Reagent (B) Antiserum Volume = 5 ml	Goat serum, anti-human Lambda Sodium azide 0.95 g/l

Optional: Kappa / Lambda Calibrator – Ref. 7809 - 5X1 ml

STABILITY AND STORAGE

The reagents are stable until the expiration date indicated on the label if stored tightly closed at 2-8°C.

Once opened, the reagents are stable at least 2 months at 2-8°C protected from light and in the absence of contamination.

Do not use over expiry date. Keep bottles closed when not in use.

REAGENT PREPARATION

The Reagents are liquid, ready to use.

PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow. Use the normal precautions required in the laboratory. Dispose of waste according to local laws.

ASSAY PROCEDURE *

Sample/Calibrators/Controls: Dilute 1:3 with saline

Wavelength:	340 nm
Temperature:	37°C
Method:	Sample Blank (A)
Sample, Control or Calibrator:	5 µl
Reagent (A):	300 µl
Reagent (B):	60 µl
Incubation R(A)/R(B):	60/180 sec
Reading:	10 sec

* Validated on BT2000 plus

EXPECTED VALUES

Serum: 110 – 240 mg/dl

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of **Kappa /Lambda Controls: High** (Ref.7807- x1 ml) and **Low** (Ref. 7808 - 1x1 ml).

PERFORMANCE

Measuring Range: 0 - 450 mg/dl.

Detection Limit: 20 mg/dl.

Sensitivity: 0.0030 ABS units / concentration unit.

Precision intra-assay:

	Level 1	Level 2	Level 3
CV %	2.91	1.33	1.09

Precision inter-assay:

	Level 1	Level 2	Level 3
CV %	---	6.77	4.34

Accuracy:

Control	Assigned Value	Result
Control 1	77.7 (66.0 - 89.4) mg/dl	78.0 mg/dl
Control 2	155 (132 - 178) mg/dl	151 mg/dl

Interferences: No interference for: Hemoglobin (10 g/l), Bilirubin (20 mg/dl) and Turbidity (5%).

Correlation against a reference method: $Y = 1.049x + 3.318$ $r = 0.99$

REFERENCES

1. Lievens, M. M., J. Clin. Chem. Clin. Biochem. 27, 519-523 (1989).
2. Boege, F. et al., Lab. Med. 13, 369-374 (1989).
3. Dati, F. et al., Lab. Med. 13, 87-90 (1989).