FIBRINOGENO

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BIBLIOGRAFIA

- 1- Hirsh J., Deykin D., Poller L.: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
- 2- International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
- 3- National Committee For Clinical Laboratory Standards: Collection, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
- 4 - Douketis Jd, Lane A, Milne J, Ginsberg Js. Trombosis Research, 1998; 92: 11-17.

FIBRINOGEN

Method based on Thrombin clotting time

REF. 1004/L 2x5ml + 2x50ml (without Calibrator and controls)

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE

Method for the quantitative determination of fibrinogen in citrated plasma.

PRINCIPLE

The thrombin clotting time fibrinogen assay is based on the method originally described by Clauss. In the presence of high concentrations of thrombin, the time required for clot formation in dilute plasma is inversely proportional to the fibrinogen concentration.

SAMPLE

Plasma in trisodium citrate 3.2 % (0.109 M).

Avoid hemolysis and contamination by tissue fluids. Centrifuge blood for 15 minutes at 1500 x g. Test within 2 hours if samples are held at 22-24°C.

KIT COMPONENTS

Bovine Thrombin (2 x 5 ml)	Buffered Bovine Thrombin	Thrombin 100 NIH Units/ml Albumin > 1% Preservative > 1%
Imidazole Buffer (2 x 50 ml)	Imidazole buffer in saline Sodium azide 0.2 g/l	Imidazole > 10 mmol/L Preservative < 1%

Optional:

Fibrinogen Calibrator REF. 1006 (1x1 ml)	Human plasma collected with sodium citrate anticoagulant
Control Plasma N REF. 1007 (1x1 ml)	Human plasma collected with sodium citrate anticoagulant < 0.4 %
Control Plasma P REF. 1008 (1x1 ml)	Human plasma collected with sodium citrate anticoagulant < 0.4 %

The reagents are stable until the expiration date indicated on the label if stored tightly closed at 2-8°C.

Erratic values, product color variations, or lack of vacuum in the vial of Bovine Thrombin, Calibrator or Controls could indicate product deterioration. However, poor control performance could also be due to other factors within the test system.

REAGENT PREPARATION

The reagents are ready to use.

Reconstitute the Fibrinogen Calibrator, Control Plasma N, and Control Plasma P with 1.0 mL of distilled water. Gently swirl and let stand for 15 minutes at room temperature. Do not invert or shake the vial vigorously. Once properly reconstituted, the calibrator and controls remain stable for 2 and 8 hours respectively at 2-8°C.

PRECAUTIONS AND WARNINGS

Use standard laboratory precautions.

The reagent may contain non-reactive components and various types of preservatives. As a precaution, avoid skin contact and ingestion. Use standard laboratory safety procedures.

Dispose of waste according to local regulations.

TEST PROCEDURE

Prepare various dilutions (1/5, 1/10, 1/15, 1/20, 1/40) of the calibrator as shown in Table 1.

Table 1. Calibration curve:

	Fibrinogen Calibrator	Imidazole Buffer	Total Volume
1/5	100 µl	400 µl	500 µl
1/10	100 µl	900 µl	1000 µl
1/15	100 µl	1400 µl	1500 µl
1/20	100 µl	1900 µl	2000 µl
1/40	100 µl	3900 µl	4000 µl

Dilute the samples and controls 1:10 with Imidazole Buffer (50 µL + 450 µL).

Treat calibrators and controls in the same way as the samples.

Prediluted samples 200 µl

Incubate 4-6 minutes at 37°C. Add:

Bovine Thrombin 100 µl

Time clot formation. Do not prewarm thrombin.

The frequency of curve preparation is partially determined by the method of clot detection used. Always prepare a new curve with each change in reagent lots, instrumentation, or when controls no longer fall within established ranges.

RESULTS INTERPRETATION

Calculate the mean of duplicate clotting times and construct a log-log curve that plots fibrinogen concentration vs. clotting time.

Plasma diluted 1:10 represents 100 % of the assigned value. The dilution factor indicates the relationship between the 1:10 dilution and other dilutions.

Example: Fibrinogen Calibrator = **286 mg/dl** (each laboratory must prepare curves with their reagents and instrumentation).

Dilution	Dilution factor	Fibrinogen (mg/dl)
1:5	10/5 = 2	286x2 = 572
1:10	10/10 = 1	286x1 = 286
1:15	10/15 = 0.67	286x0.67 = 191.6
1:20	10/20 = 0.5	286x0.5 = 143
1:40	10/40 = 0.25	286x0.25 = 71.5

Find the clotting time of controls and patient samples on the curve and read the corresponding fibrinogen value.

LIMITATIONS

- Blood must be immediately added to trisodium citrate anticoagulant and gently mixed. **EDTA and Heparin are unsuitable anticoagulants.**
- Hemolysis can cause clotting factor activation and end point detection interference. Icteric and lipemic specimens may also be inappropriate for end point detection methods.
- The sample should only contact nonwettable surfaces.
- The ratio of blood to anticoagulant is usually 9:1 and results in a citrate concentration of 10.9 to 12.9 mmol/L. This concentration must be adjusted for patients with hematocrits above 55 %.
- Freezing and thawing of plasma that contains residual cells will generate damaged cell membranes that can affect results.
- Acute inflammatory reactions can elevate circulating Factor I (Fibrinogen).
- High Fibrinogen Degradation Products (FDP) may prolong clotting times, especially when the fibrinogen level is below 150 mg/dl.
- In patients with qualitative fibrinogen abnormalities, the thrombin clotting time assay may indicate decreased fibrinogen. The quantitative fibrinogen results may be normal on these same samples if tested by other methods.
- Heparin does not interfere at therapeutic levels. However, very high heparin levels may cause low fibrinogen results. Batroxobin enzyme can be substituted for thrombin in this assay if heparin interference is suspected.
- High paraprotein levels, thrombin antibodies, and drugs that activate the fibrinolytic system can interfere with fibrinogen assays.
- The kit Fibrinogen and individual components are designed to work at 37°C. Ensure that all heating elements are functioning properly.

EXPECTED VALUES

Generally the normal reference interval is:

150 – 350 mg/dl (1.5 – 3.5 g/l)

These values should only be used as a guideline. Each laboratory should establish own Reference Range.

PERFORMANCE

Precision: Twenty consecutive determinations were performed over three consecutive days using three samples at different concentrations. A percentage CV of 1.37% (low), 0.96% (normal), and 0.79% (high) was observed

Correlation with reference method: $Y = 1.0082x - 3.4767$ $r = 0.9992$

REFERENCES

- 1- Hirsh J., Deykin D., Poller L.: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
- 2- International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
- 3- National Committee For Clinical Laboratory Standards: Collection, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
- 4 - Douketis Jd, Lane A, Milne J, Ginsberg Js. Trombosis Research, 1998; 92: 11-17