COLESTEROLO HDL

Reagente precipitante PEG 6000

REF. 0056 2x100 ml

IVD

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REF. 0056 2x100 ml

IVD

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE

Reagent for the precipitation of β -lipoproteins VLDL (very low-density lipoproteins) and LDL (low-density lipoproteins), to be used for the quantitative determination of HDL Cholesterol in serum.

PRINCIPLE

The β -lipoproteins VLDL and LDL are precipitated with a polyanionic reagent (PEG 6000) and on the supernatant, after centrifugation, the alpha fraction, HDL cholesterol (high density lipoproteins) is determined with colorimetric enzymatic method.

SAMPLE

Fresh unhemolyzed serum

HDL in the sample is stable 3 days at 2-8°C and one month at -20°C.

KIT COMPONENTS

Reagent (A) HDL Volume = 2x100 ml	Good Buffer PEG 6000 Sodium Azide	100 mmol/l 100 mmol/l 14 mmol/l
Standard CHOL Volume = 1x5 ml	Cholesterol Sodium azide	Value on label 14 mmol/l

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C and protected from light. Do not freeze. Once opened reagents are stable for 2 months at 2-8°C if contamination is avoided.

Keep bottles closed when not in use.

REAGENTS PREPARATION

Liquid Reagents, bring to room temperature (15-25°C) before use.

PRECAUTIONS AND WARNINGS

Reagents may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength:	510 nm (500 – 520)
Lightpath:	1 cm
Temperature:	37°C
Reading:	against blank reagent
Method:	Increasing End Point

PHASE 1:

Pipette:

Reagent (A) HDL	500 μ l
Sample	500 μ l

Mix on vortex and incubate for 5 minutes at room temperature, centrifuge at 2000/3000 r.p.m. for 10 minutes.

N.B. In Phase 1, the sample is diluted 1:2; therefore, it will be necessary to multiply the result obtained in Phase 2 by 2 (dilution factor).

Incomplete precipitation of VLDL and LDL will rarely occur. However, if it becomes evident that precipitation is not complete (lipid supernatant), it is recommended to increase the dilution factor of the sample and use 3 parts of HDL Reagent (A) and 1 part of sample. Then multiply the result obtained at the end of Phase 2 by 4 (dilution factor).

PHASE 2:

Determine the HDL Cholesterol in the supernatant obtained during Phase 1 using the **CHOLESTEROL SL** kit, enzymatic colorimetric method (REF. 0045 – 0035 – 0093 – 0045L), as follows:

Pipette:	Blank	Sample	Standard
Reagent (A) CHOLESTEROL SL	1000 μ l	1000 μ l	1000 μ l
Water	25 μ l		
Sample		25 μ l	
Standard			25 μ l

Mix, incubate at 37°C for 5 minutes, read the absorbance of the sample (Ax) and the standard (As).

Volumes can be proportionally modified.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

RESULTS CALCULATION

HDL (mg/dl) = $A_x/A_s \times (\text{Standard Value}) \times 2$ (dilution factor)

For the calculation of LDL Cholesterol use the following formula:

LDL (mg/dl) = Total Cholesterol – HDL – (Triglycerides/5)

This formula will be valid if the value of triglycerides is less than 400 mg/dl and the patient does not have a type III hyperlipoproteinemia.

Conversion factor: mg/dl x 0.02586 = mmol/l

EXPECTED VALUES

HDL:

Women: 42 - 88 mg/dl (1.09 - 2.28 mmol/l)

Men: 35 - 80 mg/dl (0.90 - 2.07 mmol/l)

LDL:

Adults: 66 - 178 mg/dl (1.71 - 4.60 mmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of the control (REF.6011).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 2 mg/dl.

Linearity: the method is linear up to 700 mg/dl. For higher values, dilute the sample 1:5 and multiply the result by 5.

Precision intra-assay:

	Level 1	Level 2
Mean (mg/dl)	32.1	88.9
SD	0.18	0.61
CV %	0.55	0.68

Precision inter-assay:

	Level 1	Level 2
Mean (mg/dl)	32.8	100.1
SD	0.4	1.1
CV %	1.3	1.1

Interferences: bilirubin does not interfere up to 20 mg/dl. Hemoglobin does not interfere up to 500 mg/dl. Triglycerides do not interfere up to 1000 mg/dl.

Correlation against a reference method: $Y = 0.96x + 2.5$ $r = 0.998$

REFERENCES

- Trinder P., Ann. Clin. Biochem. 6, 24 (1969).
- Kaplan LA, Pesce AJ, Clinical Chemistry, Mosby Ed. 1989.
- Vassault A. et al., Ann. Biol. Clin., 44, 686 (1986).
- Young D.S., et al., Clin. Chem. 21:1D (1975).
- Friedewald W.F., et al., Clin. Chem. 18: 449-502, 1972.