PROTEINE TOTALI

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Colorimetric Method - Biuret
Liquid Reagent ready to use

REF. 0046/2 2x 100 ml

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE
DIAGNOSTICS®

Giesse Diagnostics srl
V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18,300 - 00012 Guidonia Montecelio (RM) - Italia
Tel. +39 0774 051100 - Fax +39 0774 051111
e-mail: info@giessediagnostics.com - web site: www.giessediagnostics.com

004607
Ed. 2025/11 rev. 08

REF. 0046 4x 100 ml
REF. 0042 4x 250 ml

IVD

INTENDED USE

Quantitative determination of total proteins in serum and plasma.

PRINCIPLE

In an alkaline medium, proteins form a blue-violet complex with Cu (II) ions. The color intensity is proportional to the concentration of total proteins in the sample.

SAMPLE

Serum, plasma with heparin. Do not use hemolyzed samples.

Proteins in the sample are stable one week at room temperature (15-25 °C) and one month at 2-8°C.

KIT COMPONENTS

Reagent (A) PT Volume = 100/250	Sodium hydroxide	190 mmol/l
	Copper sulphate	8 mmol/l
	Potassium iodide	33 mmol/l
	Potassium sodium tartrate	50 mmol/l

Optional: Standard for Total Proteins REF. 8801 o 8801/6

The Standard is not included in the kits.

Standard (B) PT Volume = 5 ml	Total Proteins (Albumin) Sodium azide	Value on label
----------------------------------	--	----------------

The reagents are stable until the expiration date if stored at temperatures indicated on the label and protected from light. Once opened, the reagent are stable 2 months if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagent ready to use.

PRECAUTIONS AND WARNINGS

EC Regulation 1272/2008 (CLP):

Warning

H319 Causes serious eye irritation.

H315 Causes skin irritation.

H412 Harmful to aquatic life with long lasting effects..

EUH208 Contains: POTASSIUM IODIDE
May produce an allergic reaction.

P280 Wear protective gloves / eye protection / face protection.

P264 Wash with water thoroughly after handling.

P273 Avoid release to the environment.

P337+P313 If eye irritation persists: Get medical advice / attention.

Reagents may contain some non-reactive and preservative components.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength: 546 nm (530 – 550)

Lightpath: 1 cm

Temperature: 25, 30, 37°C

Reading: against blank reagent

Method: Increasing End Point

Sample/Reagent: 1/100

pipette:	blank	Sample	standard
Reagent (A)	1000 µl	1000 µl	1000 µl
water	10 µl		
sample		10 µl	
standard			10 µl

Mix, incubate at 25, 30, 37°C for 10 minutes and read against blank reagent the absorbance of the sample (Ax) and the standard (As).

Reaction volumes can be proportionally varied.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

RESULTS CALCULATION

Total Proteins g/dl = $Ax/As \times 6$ (standard Value)

EXPECTED VALUES

Serum/plasma: 6.6 – 8.3 g/dl

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 0.1 g/dl.

Linearity: the method is linear up to 13 g/dl. For higher values, dilute with saline 1:2 the sample and multiply the result by 2.

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	Level 1	Level 2	Level 3
Mean (g/dl)	2.58	5.08	9.75
SD	0.043	0.044	0.110
CV %	1.65	0.87	1.13

Precisione tra le serie:

	Level 1	Level 2	Level 3
Mean (g/dl)	2.65	5.36	10.12
SD	0.028	0.052	0.156
CV %	1.04	0.97	1.54

Interferences: bilirubin does not interfere up to 30 mg/dl. A triglycerides concentration exceed 300 mg/dl increases the reading.

Hemoglobin in concentration > 0.3 mg/dl increases the reading.

Correlation against a reference method:

$Y = 1.0068x - 0.0670$ $r = 0.9974$

REFERENCES

1. Vassault, A. et al. Ann. Biol. Clin., 44,686 (1986).
2. Young, D.S., et al., Clin. Chem. 21:1D (1975).
3. Kaplan LA, Pesce AJ, Clinical Chemistry, Mosby Ed. 1989
4. Kingsley G.R., Biol. Chem. 131, 197-200 (1939)

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18,300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

004607
Ed. 2025/11 rev. 08