UREA 2

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UREA 2

Modified Berthelot method (enzymatic colorimetric)

REF. 4063
2x250 ml

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by **IMQ**

INTENDED USE
Quantitative determination of urea in serum, plasma and urine.

PRINCIPLE
Urease catalyzes the hydrolysis of urea into carbon dioxide and ammonia which, in the presence of sodium nitroprusside, reacts with sodium hypochlorite and sodium salicylate, forming a green compound whose intensity is proportional to the concentration of urea present in the sample.

SAMPLE
Serum, plasma, urine.
Avoid anticoagulants containing ammonium salts and fluoride.
The urine should be diluted 1:100 with saline.
Urea in the sample is stable 3 days at 2-8 °C.

KIT COMPONENTS		
Reagent (A) Volume = 240 ml	Good Buffer Sodium Nitroprusside Sodium Salicylate	50 mmol/l 3.5 mmol/l 90 mmol/l
Reagent (B) Volume = 10 ml	Good Buffer Urease Sodium azide	20 mmol/l 36000 U/l 3 mmol/l
Reagent (C) Volume = 19 ml	Sodium hydroxide Sodium hypochlorite	200 mmol/l 500 mmol/l
Standard Volume = 5 ml	Urea Sodium azide	Value on label

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C. Once opened reagents are stable for 2 months at 2-8°C if contamination is avoided.
Keep bottles closed when not in use.

REAGENTS PREPARATION
Add the content of a **Reagent (B)** bottle in a **Reagent (A)** bottle. Mix gently.
The working solution (A+B) is stable 5 days at 15-25°C and 3 weeks at 2-8°C.
Pour in the appropriate plastic bottle the **Reagent (C)** and add 230 ml of distilled water. This working solution is stable 3 months at 2-8°C.

PRECAUTIONS AND WARNINGS
EC Regulament 1272/2008 (CLP):
Reagent (A) :
EUH210 Safety data sheet available on request.
Reagent (C) :

Warning

H319 Causes serious eye irritation

H315 Causes skin irritation.

H412 Harmful to aquatic life with long lasting effects.

P280 Wear protective gloves / eye protection / face protection.

P337+P313 If eye irritation persists: get medical advice

P264 Wash with water thoroughly after handling.

P273 Avoid release to the environment.

Reagents may contain some non-reactive and preservative components.
Use the normal precautions required in the laboratory.
Dispose of waste according to local laws.

PROCEDURE	
Wavelength:	578 nm
Light path :	1 cm
Temperature:	37°C
Reading:	against blank
Method:	End Point

pipette:	Blank	Sample	Standard
Reagent (A+B)	1000 µl	1000 µl	1000 µl
Water	10 µl		
Sample		10 µl	
Standard			10 µl
Mix, leave in a water bath for 5 minutes at 37 °C, add:			
Working solution (C)	1000 µl	1000 µl	1000 µl
Mix, incubate for 5 minutes at 37°C and read against blank the absorbance of the sample (Ax) and the standard (As), within 5 minutes.			

Reaction volumes can be proportionally varied.

RESULTS CALCULATION
Serum/plasma:
Urea mg/dl = $A_x/A_s \times 50^*$

Urine
Urea mg/dl = $A_x/A_s \times 50^* \times 100$ (dilution)

24 hours Urine:
Urea g/24h = $\frac{\Delta A_x/\Delta A_s \times 50^* \times 100 \times \text{Urine Volume (dl)}}{1000}$

*Standard Value on label

Conversion factor: mg/dl x 0.1665 = mmol/l

EXPECTED VALUES			
Serum/plasma:	10 – 50	mg/dl	(1.7 – 8.3 mmol/l)
Urine:	847 – 2967	mg/dl	(141 – 494 mmol/l)
24h Urine:	10 – 35	g/24h	(170 – 580 mmol/24h)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL
You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE
Sensitivity: the sensitivity of the method: 1.5 mg/dl.
Linearity: the method is linear up to 300 mg/dl. For higher values, dilute with saline 1:10 the sample and multiply the result by 10.
Precision intra-assay:

	Level 1	Level 2
Mean (mg/dl)	19.72	171.91
SD	1.12	1.21
CV %	2.98	3.81

Precision inter-assay:		
	Level 1	Level 2
Mean (mg/dl)	18.71	175.68
SD	1.33	4.00
CV %	5.25	2.16

Interferences: bilirubin does not interfere up to 30 mg/dl. Hemoglobin does not interfere up to 100 mg/dl.

Correlation against a reference method: $Y = 0.986x - 0.64$ $r = 0.985$

REFERENCES

- Krieg M. et al., J Clin Chem Clin Biochem 24:863 (1986)
- Chaney, et al. Clin. Chem. and al., 8, 130 (1962).
- Vassault A. et al., Ann. Biol. Clin., 44, 686 (1986).
- Mackay EM, Mackay LL. J Clin Invest 4:295 (1927)