CK-MB

Metodo cinetico UV

REF. 0083

40+10 ml

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BIBLIOGRAFIA

1. Abbot B et al. Creatinine Kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116
2. Gerhardt W et al. Creatine Kinase B-Subunit activity in serum after immunohinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young, D.S., Effects of drugs on Clinical Lab. Tests, 4th Ed., AACC Press, 1995.
4. Young, D.S., Effects of disease on Clinical Lab. Tests, 4th Ed., AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd Ed. AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd Ed. AACC 1995.
7. Mathieu M. et coll. Recommendation pour la mesure de la concentration catalytique de la creatinine kinase dans la serum humain. Ann. Biol. Clin., 40, (1482), 87.

CK-MB Kinetic UV Method		 
REF. 0083	40+10 ml	

COMPANY WITH
QUALITY SYSTEM
ISO 9001 - ISO 13485
Certified by IMQ

GIESSE[®]
DIAGNOSTICS

INTENDED USE
Quantitative determination of CK-MB in serum.

PRINCIPLE
CK-MB consists of the subunits CK-M and CK-B.
Specific antibodies against CK-M inhibits completely CK-MM activity (main part of the total CK activity) and CK-M subunit of CK-MB. Therefore only CK-B activity is measured, which is half of the CK-MB activity.
The activity of CK-MB is obtained multiplying the activity of CK-B by 2

SAMPLE
Serum, plasma with heparin. Do not use hemolyzed samples.
CK activity in serum is unstable and is rapidly lost during storage. CK is inactivated both by bright daylight and by increasing specimen pH owing to loss of carbon dioxide; accordingly, specimens should be stored in the dark in tightly closed tubes. The CK is subject to thermal denaturation, rapidly cool the sample at 4°C after collection.

KIT COMPONENTS			
Reagent (A) CK-MB Volume = 40 ml	Imidazole buffer, pH 6.7	100 mmol/l	
	D-Glucose	20 mmol/l	
	N-Acetyl-L-cysteine	20 mmol/l	
	Magnesium acetate	10 mmol/l	
	NADP	2.0 mmol/l	
	EDTA	2.0 mmol/l	
	Hexokinase	> 2.5 U/ml	
	Anti-CK-M antibody with inhibitory capacity	8000 U/l	
Reagent (B) CK-MB Volume = 10 ml	ADP	2.0 mmol/l	
	AMP	5.0 mmol/l	
	Creatine phosphate	30 mmol/l	
	Diadenosine pentaphosphate	10 µmol/l	
	G-6-PDH	> 1.5 U/ml	
	Preservatives		

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C and protected from light.

Do not use over expiry data.

Once opened reagents are stable for 2 months at 2-8°C if contamination is avoided. Keep bottles closed when not in use.

REAGENT PREPARATION
Add 1 Volume of Reagent (B) to 4 Volumes of Reagent (A). The working solution (A+B) is stable 20 days at 2-8°C.

Gently mix before use.

PRECAUTIONS AND WARNINGS
Reagent (A) – Contains Acetic Acid (CAS 64-19-7)

Signalwords: **Warning**

H360D May damage the unborn child.

P201 Obtain special instructions before use.

P280 Wear protective gloves/ protective clothing / eye protection / face protection.

P308+P313 IF exposed or concerned:

Get medical advice / attention.

P202 Do not handle until all safety precautions have been read and understood.

P405 Store locked up.

P501 Dispose of contents / container to an approved waste disposal plant.

The reagent may contain non-reactive components and preservatives of various kinds. It is recommended to handle with care, avoid contact with skin and ingestion. Use standard laboratory safety precautions. Dispose of waste in accordance with local regulations.

PROCEDURE	
Wavelength:	340 nm
Lightpath :	1 cm
Temperature:	37 °C (25 °C o 30 °C)
Reading:	against distilled water
Method:	increasing kinetic

Reaction volumes can be proportionally varied.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

pipette:	sample	blank	calibrator
Working solution (A+B)	1000 µl	1000 µl	1000 µl
sample	40 µl	-	-
water	-	40 µl	-
calibrator	-	-	40 µl

Gently mix and incubate at 37 °C (25 °C or 30 °C) for 5 minutes, read the initial absorbance (A1) of the sample, and repeat the absorbance reading after 5 minutes (A2). Calculate the difference $\Delta A = A1 - A2$. The same procedure applies to the calibrator.

RESULTS CALCULATION

CK-MB activity (U/L) = $(\Delta A_{\text{sample}} / \Delta A_{\text{calibrator}}) \times \text{calibrator concentration}^*$

Percentage of CK-MB activity in the sample:

% CK-MB = (CK-MB / Total CK) × 100

To calculate total CK activity, use the CK-NAC SL Kits (Ref. 0018 - 4054 - 4052).

(*) The use of the CK-MB Calibrator (Ref. 6021) is recommended for this kit.

EXPECTED VALUES

CK-MB
< 24 U/l at 37 °C
< 15 U/l at 30 °C
< 10 U/l at 25 °C

% CK-MB
< 6 %

The probability of myocardial infarction is high under the following conditions:

		U/l at 25 °C	U/l at 30 °C	U/l at 37 °C
CK – MB		> 10	> 15	> 25
CK	Men	> 80	> 130	> 190
	Women	> 70	> 110	> 167
% CK – MB		6 – 25 %		

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control serum: CK-MB Control (Ref. 6020).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 1,9 U/L.

Linearity: the method is linear up to 1000 U/L. For higher values, dilute the sample 1:1 and multiply the result by 2.

Precision intra-assay:

	Level 1	Level 2
Mean (U/L)	172.1	776.4
CV %	2.83	1.73

Precision inter-assay:

	Level 1	Level 2
Mean (U/L)	165.4	740.2
CV %	3.37	2.06

Interferences: No significant interference up to 2,5 g/l of triglycerides. Glucose does not interfere up to 7g/l. Hemoglobin does not interfere up to 1,25 g/l. Bilirubin does not interfere up to 600 µmol/l.

Correlation against a reference method: $Y = 0.976x - 0.269$ $r = 0.999$

REFERENCES

1. Abbot B et al. Creatinine Kinase. Kaplan A et al. Clin Chem The C.V. Mosby Co. St Louis. Toronto. Princeton 1984: 1112-116
2. Gerhardt W et al. Creatine Kinase B-Subunit activity in serum after immunoinhibition of M-Subunit activity. Clin Chem 1979;(25/7): 1274-1280.
3. Young, D.S., Effects of drugs on Clinical Lab. Tests, 4th Ed., AACC Press, 1995.
4. Young, D.S., Effects of disease on Clinical Lab. Tests, 4th Ed., AACC 2001.
5. Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd Ed. AACC 1999.
6. Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd Ed. AACC 1995.
7. Mathieu M. et coll. Recommandation pour la mesure de la concentration catalytique de la creatinine kinase dans la serum humain. Ann. Biol. Clin., 40, (1482),87.

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 05111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

008307

Ed. 2025/11 rev.04